A stylized illustration of space. In the upper left, a large, grey, cratered moon is visible. In the upper right, a reddish-brown planet (Mars) is shown. A bright yellow star is in the lower left, and a pink star is in the top right. A blue comet with a long tail is streaking across the sky. In the bottom left, a portion of a space shuttle is visible, with its orange external tank and white boosters. The title text is centered in the upper right area.

VERIFICATION OF NASA GALACTIC COSMIC RAY SIMULATOR FOR LARGE ANIMAL MODELS

Shirin Rahmanian¹, Tony C. Slaba²

¹*National Institute of Aerospace, Hampton, VA 23666*

²*NASA Langley Research Center, Hampton, VA 23681*



Outline

- Galactic Cosmic Ray Simulator (GCRsim) Beam Overview
- Transport Studies
- Digital Phantoms
- Simulation Results of GCRsim Beam in the Three Different Phantoms
- Summary and Conclusions



GCRsim Beam Overview: Background

- Space radiation poses multiple important health risks for astronauts
 - Cancer
 - Cardiovascular disease
 - Damage to Central Nervous System
- For long duration mission beyond low Earth orbit (LEO) risks mainly arise from galactic cosmic rays (GCR)
- Ground-based experiments will help to mitigate risks and reduce uncertainties



GCRsim Beam Overview: Objective


- Previous ground-based radiation studies mainly utilized single mono-energetic beams
 - Collectively improved our understanding of underlying biological mechanism
 - Poor analog for the complete space environment
- The GCR simulator (GCRsim) was developed at the NASA Space Radiation Laboratory (NSRL) to better represent the complex mixed field environment in space
- GCRsim is intended simulate radiation environment as seen by astronauts in deep space
 - Study health effects of GCR
 - Improve risk projections
 - Countermeasure development and testing

NASA's NSRL facility in Brookhaven, NY




GCRsim Beam Overview: Reference Field

- A single reference field was defined to approximate deep space environment
 - Female BFO (blood forming organ) behind 20 g/cm² spherical aluminum shielding during solar minimum



Reference field energy spectra for neutrons, hydrogen, and helium [Slaba et al. 2016]



Differential LET spectra of reference field with and without contributions from hydrogen and helium [Slaba et al. 2016]



GCRsim Beam Overview: Beam Definition

- The GCRsim at NSRL is designed to deliver deep space, shielded tissue environment to biological targets in a laboratory setting
- 33 mono-energetic beams of varying energies with ion species consisting of H, He, C, O, Fe, Si and Ti
- Sequential beam delivery reproducing the space environment over the full range of LET

Ion	Energy (MeV/n)	Range (cm)	LET (keV/μm)	Dose (mGy)
¹ H	100	<i>Polyethylene degrader to</i>		
¹ H	150	15.9	0.54	35.0
¹ H	250	38.1	0.39	68.9
¹ H	1000	326.6	0.22	123.6
⁴ He	100	<i>Polyethylene degrader to</i>		
⁴ He	150	16.0	2.17	7.5
⁴ He	250	38.3	1.56	16.4
⁴ He	1000	327.8	0.88	24.9
¹² C	1000	110.1	7.95	11.7
¹⁶ O	350	17.0	20.8	15.4
²⁸ Si	600	22.7	50.2	8.1
⁴⁸ Ti	1000	32.5	109.5	4.5
⁵⁶ Fe	600	13.1	175.1	4.1
Total				500.0

GCRsim beam definition at NASA
[Simonsen et al. 2020]

Ion	Energy (MeV/n)	Range (cm)	LET (keV/μm)	Dose (mGy)
¹ H	20.0	0.43	2.59	30.4
¹ H	23.3	0.56	2.29	6.7
¹ H	27.2	0.75	2.02	7.4
¹ H	31.7	0.98	1.79	8.0
¹ H	37.0	1.30	1.58	8.7
¹ H	43.2	1.72	1.39	9.3
¹ H	50.3	2.26	1.23	10.0
¹ H	58.7	2.99	1.09	10.6
¹ H	68.5	3.95	0.97	11.1
¹ H	79.9	5.20	0.86	11.2
¹ H	100.0	7.76	0.73	27.2

Ion	Energy (MeV/n)	Range (cm)	LET (keV/μm)	Dose (mGy)
⁴ He	20.0	0.43	10.34	11.0
⁴ He	23.3	0.57	9.14	2.1
⁴ He	27.2	0.75	8.06	2.2
⁴ He	31.7	0.99	7.12	2.3
⁴ He	37.0	1.31	6.29	2.5
⁴ He	43.2	1.73	5.56	2.6
⁴ He	50.3	2.28	4.92	2.7
⁴ He	58.7	3.01	4.36	2.7
⁴ He	68.5	3.97	3.86	2.7
⁴ He	79.9	5.23	3.43	2.7
⁴ He	100.0	7.81	2.90	6.1




Transport Studies

- The GCRsim is suitable for animal models such as mice and rats to represent internal radiation environment seen at critical organ locations within the human body
- Radiation transport studies using phantom models of mouse (Digimouse) and rat (Digirat) have been previously completed [Simonsen et al. 2020]
 - Geant4 Monte Carlo code was used to simulate the GCRsim in the digital phantoms
- Verified key physical parameters of the GCRsim beams
 - Homogeneous internal dose distribution across radiosensitive tissues
 - Reproduces the dose and fluence spectra of the reference field as function of linear transfer energy (LET)



Digital Phantoms


- Digital phantoms are voxelized 3D models of an animal created from CT or MR images
 - Digimouse: Digital model of 28 g male mouse
 - Digirat: scaled model of Digimouse resulting in a 754 g rat
 - Digipig: Digital model of a 35 kg male minipig
- Each voxel identifies important radiosensitive tissues
 - Quantities of interest were calculated in these radiosensitive organs
- Simulations performed with each of these phantoms
 - Geant4 Monte Carlo simulation
 - 500 mGy GCRsim beam dose
 - Isotropic irradiation conditions



Digipig, Digirat and Digimouse shown side by side for comparison



Simulation Results: Dose Distribution




- Relatively homogenous dose distribution was seen throughout all three animal phantoms
- 95% of the voxel doses were within 6%, 7%, and 8% of the mean values in Digimouse, Digirat, and Digipig, respectively



Simulation Result: Tissue Dose



- No major variation in tissue doses was observed
 - All values compared well with beam dose of 500 mGy
 - Bone dose slightly lower in all animals due to density difference
- Digimouse: mean soft tissue dose was 4.6% higher than the beam dose
 - Bone dose was 6% lower
- Digirat: mean soft tissue dose 2% lower than the beam dose
 - Bone dose was 11% lower
- Digipig: mean soft tissue dose about 10% lower than the beam dose
 - Bone dose was 20% lower



Side by side comparison of simulated dose in major radiosensitive organs of Digimouse, Digirat and Digipig.



Simulation Result: Fluence and Dose Spectra



Comparison of the cumulative fluence (on the left), and cumulative dose (on the right) as a function of LET within radiosensitive organs of the three different animal phantoms

- Simulated spectra in all three animal models compare well with reference field
- The difference between the phantoms and the reference field was less than 15% in LET domain that contributes most heavily to dose



Summary and Conclusion

- It was shown from the Monte Carlo simulations that the relevant radiation quantities in larger animals such as minipigs have comparable values to rodent models and the reference field.
- The differences observed were within expectations due to the increased mass and most likely not substantial enough to significantly change biological responses.
- These results suggest that no major modifications may be required for the existing GCRsim beam to support studies with large animal models.
- Nevertheless, further analysis with large animals may be needed to evaluate the radiation field in the internal organs for specific experimental design considerations.



References

- [Simonsen et al. 2020] Simonsen, L. C., T. C. Slaba, P. Guida, and A. Rusek. "NASA's first ground-based galactic cosmic ray simulator: Enabling a new era in space radiobiology research." PLoS Biology 18, no. 5 (2020): e3000669.
- [Slaba et al. 2016] Slaba, T. C., S. R. Blattnig, J. W. Norbury, A. Rusek, and C. La Tessa. "Reference field specification and preliminary beam selection strategy for accelerator-based GCR simulation." Life Sciences in Space Research 8 (2016): 52-67

